

Mobi-walk

Beach and all-terrain rolling walker



Ideal for walking on various types of soft ground (dry and wet sand, grass, earth, gravel, etc.), the Mobi-walk walker enables people with limited abilities to move around independently. Its two rear wheels are equipped with brakes to ensure user safety.



ADVANTAGES



- > Perfect wheels for soft grounds
- > Safety and comfort: brakes on back wheels, padded and foldable backrest
- > The pleasure of sitting down to enjoy nature
- > Foldable: easily transportable
- > Adjustable height: 5 positions

Technical specifications

Colors	Silver frame, grey wheels and black handles			
Materials	Mobi-walk frame is made of anodized aluminum and handles are made of foam.			
Dimensions	Length	Width	Width between handles	Weight & weight capacity
	91,2 cm 36"	90 cm 35.4"	50,5 cm 20"	10,2 kg 22.5 lbs 125 kg 275.5 lbs
	Seat height	Handle width	Seat dimensions	Wheel diameter
	55 cm 22.6"	86cm - 95cm du sol 34" x 37.4" on the ground	41cm x 28cm 16"x11"	Arrière : 24cm x 12.3cm Avant : 30cm x 18cm Back : 9.4" x 4.8" Front : 11.8" x 7"

General Characteristics

The frame of the Mobi-walk walker is made of aluminum with stainless steel screws. The wheels with bearings are unbreakable. The handles are ergonomic and adjustable, and there is a brake for each wheel. The walker is foldable and can support a weight of 125 kg | 275.5 lbs.

Ideal for walks on several types of grounds (dry and wet sand, grass, soil, gravel, etc), Mobi-walk rollator allows people with reduced mobility to move independently.

Its two rear wheels are equipped with brakes which ensure safety for the user. Mobi-walk all-terrain rollator adjusts to each person and allows you to immerse yourself in nature independently. This will make it easier to access various soft grounds.

Easy installation & low maintenance

QUICK INSTALLATION

The setup of the Mobi-walk is simple and quick. First, fully unfold the walker so that it is stable. Next, adjust the handles to the appropriate height to fit the user and ensure a comfortable and secure posture.

Once the adjustments are made, it is important to check that the brakes are working properly and that all fastenings are securely tightened. Before use, the seat can be installed if needed. These checks ensure safe and comfortable use of the Mobi-walk on all types of terrain.

EASY MAINTENANCE

Maintenance of the Mobi-walk is limited to regular cleaning after use, especially after exposure to sand, with drying before storage. It is important to check before each use that the wheels and brakes are functioning properly and that all fastenings are securely tightened.

Periodic inspection of parts for wear and storing the walker in a dry place help ensure the safety and longevity of the equipment.



What are our additional products?

Mobi-mat® Accessibility mat



Mobi-roller™ Mats dispenser



Mobi-com™ Advertising mat



Sunbed With adapted seat height

